

Welch Allyn Green Series™ 777 Wall Transformer

FOR VETERINARIANS



(instrument heads
not included)

Increased energy savings to help you improve your bottom line.*

- Improve exam efficiency with no on/off switch and save staff time without the need to turn unit(s) on/off each day
- Allows for reliable patient diagnosis in a standardized solution
- Meet your exam needs with handle-based rheostats for quick light-intensity adjustments
- Reliable optical sensor automatically turns handles on when they are removed from the cradles
- Easy to maneuver, contoured handles prepare you for any type of exam
- Coiled cords (12' extended) are guaranteed for 10 years
- Quickly and easily upgrade from the Welch Allyn 767 Wall Transformer

**More energy efficient than the Welch Allyn 767 Wall Transformer*

Specifications

Product Description:	Welch Allyn Green Series™ 777 Wall Transformer
Product Type:	Wall Transformer
Unit Dimensions:	8.07" x 9.25" x 4.92" (205 mm x 235 mm x 125 mm)
Weight:	2.75 Lbs (1.25 Kg)
Power Source:	100–240 V, 50/60 Hz, 0.18 A max 5 V, 1 A
Warranty:	Five-year limited on the Wall Transformer, 10 years for the Handles and Cords

Ordering Information

- #77710 Green Series 777 Wall Transformer with 8' Twin USB Power Cable
#104714 18" Twin USB Power Cable

www.welchallyn.com/dvm



WelchAllyn®

Advancing Frontline Care™

Welch Allyn Diagnostic Instruments

Configure the Welch Allyn 777 Wall Transformer to meet your practice's needs.

Welch Allyn Oscopes & Ophthalmoscopes

- Standardize or choose the technology you need for each facility
- Upgrade to the Welch Allyn LED Lamp for greater energy efficiency and lower cost of ownership



#23862

Welch Allyn MacroView™ Otoloscope

- Almost twice the field of view and 30% greater magnification than traditional otoscopes
- Allows for simplified diagnosis of Otitis Media through enhanced visualization of the tympanic membrane



#11800-V and #11820

Welch Allyn PanOptic™ Ophthalmoscope

- Fast, easy entry into undilated pupils with 5X larger view of the fundus for a streamlined exam
- 25° field of view vs. the standard 5° field of view provides optimal visualization



#20260

Welch Allyn Pneumatic Otoloscope

- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Easier instrumentation to allow a faster and more comfortable exam for the patient



#11720

Welch Allyn Coaxial Ophthalmoscope

- 18 unique aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters



#21760

Welch Allyn Operating Otoloscope

- Open system for convenience during procedures
- Rotatable lens and speculum for ease of use

Accessories

#23862	MacroView Veterinary Otoloscope with Specula
#21760	3.5 V Operating Otoloscope with Specula
#20260	3.5 V Pneumatic Otoloscope with Specula
#11800-V	PanOptic Ophthalmoscope with 3 spot sizes: micro, small and large; with add-on magnifying lens
#11820	PanOptic Ophthalmoscope with 3 spot sizes: micro, small and large; split aperture, red-free and cobalt blue filters; with add-on magnifying lens
#11710	3.5 V Standard Ophthalmoscope
#11720	3.5 V Coaxial Ophthalmoscope
#27000	Larynx Illuminator, complete
#27050	Nasopharynx Illuminator, complete
#41100	Finnoff Ocular Transilluminator
#43300	Curved All-purpose Transilluminator
#26038	Bivalve Specula

For more information or to order, contact your local Welch Allyn veterinary distributor, visit welchallyn.com/dvm or call 1.800.535.6663



WelchAllyn®

Advancing Frontline Care™